

Efficiency

The GM Pancake Probe efficiency is shown below. In a recent performance check, the numbers shown represent typical results obtained:

¹⁴ C	5%
⁹⁹ Tc	12%
¹³⁷ Cs	24%
⁹⁰ Sr	59%
³⁶ Cl	26%
²⁴¹ Am	8%
¹²⁹ I	2%
²³⁰ Th	15%
²³⁹ Pu	12%

NOTE: The efficiency formula used to calculate the % Efficiency is: $Eff.\% = (CPM \times 100) / DPM$